



## QUOTATION

Braun Industries, Inc  
 1170 Production Drive  
 Van Wert, OH 45891  
 800-22BRAUN

F650: 182" WB  
 Overall Length: 297.5"  
 Overall Width: 98"  
 Overall Height: 112"  
 Module Length: 170",

**Quote No:** 10121-0001

09/08/2011

Page 1

PART NO	DESCRIPTION	QTY
	== Super Chief 2012 Base Type I / 170" Mod - 12.000 08/26/11 ==	1
00-02-0141	TYPE I - AMBULANCE- FORD F650 - 73" INTERIOR HEADROOM CHASSIS PACKAGES	1
10-01-0182	CHASSIS - FORD F-650 CHASSIS, 182" WB 4 x 2	1
20-00-0018	OPTIONAL CHASSIS ITEMS - F650	1
20-00-0101	FUEL TANK - FULL UPON FACTORY RELEASE	1
20-10-0101	AIR RIDE SYSTEM, AUTOMATIC LOWERING FEATURE - OEM	1
20-10-0115	AIR DRYER	1
20-10-0214	BATTERY - THREE SYSTEM - F650	1
20-10-0402	CONDENSER AUXILIARY - BOTTOM MOUNT	1
20-10-0404	SECONDARY COMPRESSOR	1
20-10-0603	ENGINE BLOCK HEATER - OEM PLUG	1
20-20-0127	FLOOR CONSOLE, ANGLED FRONT - REMOVABLE PLATES - ALUMINUM	1
20-20-0251	HEAT SHIELDS	1
20-30-0204	MIRRORS, EXTERIOR OEM POWER TELESCOPING W/REMOTE & HEAT	1
20-30-0801	EXTERIOR CAB STEPS - NFPA DIAMOND PLATE	1
20-40-0206	WHEEL SIMULATORS, STAINLESS STEEL WITH VALVE STEM EXTENSIONS	1
30-00-0002	OPTIONAL ELECTRICAL ITEMS	1
30-06-0002	ELECTRICAL SYSTEM - WELDON V-MUX, MICRO PROCESSOR BASED	1
30-06-0026	WARNING DISPLAY - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0050	DAYTIME RUNNING LAMPS - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0076	MASTER CONTROL CONSOLE (MCC) - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0101	ATTENDANT CONTROL PANEL (ACP) - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0126	ELECTRICAL COMPARTMENT (PDQ) - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0150	SPARE CIRCUIT - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0175	WIRING HARNESS CONNECTORS - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0200	RADIO TIE-IN POINTS - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0225	EMERGENCY MASTER SWITCH - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0250	CONTROL PANEL SWITCHES - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0275	MODULE POWER SWITCH - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0285	MODULE DISCONNECT - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0300	VOLTMETER - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0325	FLASHER - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0350	FLASHING HEADLIGHTS - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0375	PARK OVERRIDE - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0400	LOAD MANAGER-SQUENTIAL SWITCHING SYSTEM - MULTIPLEX ELECTRIC	1
30-06-0425	WARNING LIGHTS SWITCH - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0450	FRONT CLEAR DISABLE SWITCH - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0475	CURBSIDE SCENE LIGHTS "ON" WITH CURBSIDE DOOR - MULTIPLEX EL	1
30-06-0500	REAR SCENE LIGHTS "ON" WITH REAR DOORS OPEN - MULTIPLEX ELEC	1
30-06-0525	DISABLE SWITCH FOR REAR SCENE LIGHTS - MULTIPLEX ELECTRICAL	1
30-06-0550	DISABLE SWITCH FOR SIDE SCENE LIGHTS - MULTIPLEX ELECTRICAL	1

PART NO	DESCRIPTION	QTY
30-06-0575	REVERSE (TRANSMISSION) FUNCTIONS - MULTIPLEX ELECTRICAL	1
30-10-0030	BACKUP ALARM	1
30-10-0050	BATTERY CONVERTER / CHARGER	1
30-10-0110	ELECTRICAL CIRCUITS	1
30-10-0190	GROUND FAULT CIRCUIT INTERRUPTER (GFCI)	1
30-10-0310	LIGHT, ELECTRICAL COMPARTMENT (PDQ)	1
30-10-0410	LIGHT OSS COMPARTMENT - INCANDESCENT	1
30-10-0472	LIGHTS, DOME HALOGEN (10) - MULTIPLEX ELECTRIC SYSTEM	1
30-10-0494	LIGHTS, FLUORESCENT (2) CURBSIDE - ANGLED	1
30-10-0496	LIGHTS, FLUORESCENT (3) STREETSIDE -ANGLED	1
30-10-0527	OUTLETS, 125VAC - THREE (3) INTERIOR DUPLEX	1
30-10-0543	OUTLETS, 12VDC - TWO (2) INTERIOR POWERPOINT	1
30-10-0657	RADIO ANTENNA BASES - TWO (2)	1
30-10-0721	RECEPTACLE SHORELINE INLET-20 AMP, NON-AUTO EJECT	1
30-10-0753	SIREN, WHELEN 295-HFSA7-REMOTE WITH MICROPHONE	1
30-10-0804	SIREN SPEAKERS, CAST PRODUCTS - IN FRONT BUMPER	1
30-10-0901	SYSTEM ACTIVE AND CHECK OUT LIGHTS - 5 MINUTES	1
30-25-0001	OPTIONAL MISC EXTERIOR LIGHTING ITEMS	1
30-25-0102	LIGHTS, ICC / LED CLEARANCE - WITH CHROME BEZELS ON MOD SIDE	1
30-25-0103	LIGHTS, ICC / LED CLEARANCE - WITH CHROME BEZELS ON FRONT	1
30-25-0104	LIGHTS, ICC / LED CLEARANCE - WITH CHROME BEZELS ON REAR	1
30-25-3107	LIGHTS, FRONT GRILLE - (2) WHELEN 700 SUPER LED - RED	1
30-25-4105	LIGHTS, FRONT MODULE WARNING -WHELEN (2)900 LED/(1)700 HALO	1
30-26-3002	LIGHTS, REAR TAIL - WHELEN 600 LED STOP, TURN, HALOGEN BK-UP	1
30-26-6010	LIGHTS, REAR MODULE WARNING-WHELEN (5) 900/ LED & HALO SCENE	1
30-26-7002	LIGHT, LICENSE PLATE LED	1
30-27-1112	LIGHTS, INTERSECTION FRONT - WHELEN 700 SUPER LED - RED	1
30-27-2007	LIGHTS, SIDE MODULE SCENE - WHELEN 900 HALOGEN - CLEAR (4)	1
30-27-3005	LIGHTS, SIDE MODULE WARNING - WHELEN 900 SUPER LED-RED	1
39-00-0001	MODULE BODY CONSTRUCTION - SEAMLESS	1
39-10-0003	CONNECTIONS, MODULE	1
40-00-0001	OPTIONAL MODULE / EXTERIOR ITEMS	1
40-00-0002	ALUMINUM PARTS - MODULE EXTERIOR	1
40-00-0010	CORROSION MINIMIZATION - MODULE	1
40-00-0020	DOOR CONSTRUCTION - GASKET PLACEMENT	1
40-00-0030	EXTERIOR COMPARTMENTS DESIGN	1
40-00-0050	EXTERIOR COMPARTMENTS, VENTING - ADDITIONAL	1
40-00-0055	EXTERIOR COMPARTMENTS, VENTING - OXYGEN STORAGE	1
40-00-0060	EXTERIOR COMPARTMENTS - SWEEP-OUT	1
40-00-0070	HANDLES, EBERHARD - CHROME FREE FLOATING TYPE	1
40-00-0080	INSULATION - POLYURETHANE SPRAY FOAM	1
40-00-0081	COATING, ACOUSTIC - AUDIOGUARD	1
40-00-0091	THRESHOLD, STAINLESS STEEL - EXTERIOR STORAGE COMPARTMENT	1
40-01-0703	OUTSIDE STORAGE #1 - STREETSIDE FORWARD - NO SHELF	1
40-02-0704	OUTSIDE STORAGE #2 - STREETSIDE CENTER - ONE (1) SHELF	1
40-03-0703	OUTSIDE STORAGE #3 - STREETSIDE REAR - (1) SHELF & TALL DOOR	1
40-04-0700	OUTSIDE STORAGE #4A BOARD/SCOOP, #4B O2 CYLINDER - RT. REAR	1
40-05-0701	OUTSIDE STORAGE #5 - CURBSIDE FORWARD	1
40-10-0066	DOOR SLIDING - MODULE CURBSIDE PASSAGE	1
40-10-0088	DOORS, REAR PASSAGE	1
40-10-0103	FENDERETTES - BRIGHT FINISH ALUMINUM	1
40-10-0120	HOLD-OPEN, REAR DOORS - CPI GRABBERS	1

PART NO	DESCRIPTION	QTY
40-10-0190	LOCKS, MANUAL (NO POWER DOOR LOCKS)	1
40-10-0200	MUD FLAPS, REAR - BLACK RUBBER	1
40-10-0261	REAR STEP/BUMPER -FLIP-UP CENTER SECTION - ANODIZED	1
40-10-0281	RUB RAILS, ANODIZED ALUMINUM, 5/8" BLACK REFLECTIVE MATERIAL	1
40-10-0304	STEP FIXED - AT CURBSIDE DOOR	1
40-10-0321	STEPWELL, CURBSIDE DOORWAY - MID STEP AREA	1
40-10-0340	STONE GUARDS, MODULE FRONT - ANODIZED DIAMOND PLATE	1
40-10-0360	THRESHOLD, REAR DOORWAY - ANODIZED DIAMOND PLATE	1
40-10-0381	TOW HOOKS, MODULE REAR	1
40-10-0505	WINDOWS, PASSAGE DOOR - (CURB SLIDING, REAR FIXED)	1
50-00-0202	ROUNDED INTERIOR CORNERS	1
50-01-0100	OPTIONAL CABINET ITEMS	1
50-01-0105	CABINET CONSTRUCTION - INTERIOR- CABINET GRADE PLYWOOD	1
50-01-0200	CABINET DOORS - BRONZE ACRYLIC	1
50-01-0402	LATCHES CABINET- BLACK	1
50-05-0109	CABINET FACING / WALL LAMINATE - ROCK OF AGE SUEDE #AW141	1
50-10-0139	CABINET - LEFT WALL, UPPER REAR #1 - FIXED FRAME	1
50-10-0212	CABINET - LEFT WALL, UPPER OVER PANEL #2 - FIXED FRAME	1
50-10-0226	LEFT WALL - ATTENDANT CONTROL PANEL - HINGED PANEL	1
50-10-0312	CABINET - LEFT WALL, UPPER FRONT#3 - FIXED FRAME	1
50-10-0412	CABINET - LEFT WALL, LOWER REAR #4 - INSIDE/OUTSIDE	1
50-10-0524	CABINET - LEFT WALL, LOWER UNDER REAR ACTION AREA #5	1
50-10-0614	PRIMARY ACTION AREA	1
50-10-0642	SECONDARY ACTION AREA	1
50-10-0901	COUNTER TOP - LAMINATED	1
50-20-0010	CABINET DOORS - LOWER FRONT WALL - ETCHED ACRYLIC	1
50-20-0276	CABINET - FRONT WALL, LOWER #1	1
50-20-0504	PARTITION DOORWAY - WITH HINGED DOOR	1
50-30-0108	CABINET - RIGHT WALL, REAR WITH (1) DOOR FACING AISLE #1	1
50-50-0001	OPTIONAL SEATING ITEMS	1
50-50-0060	SQUAD BENCH STORAGE (COMPARTMENT #2)	1
50-50-0100	ATTN SEAT - DELUXE HIGH BACK, ON SLIDERS, SWIVEL BASE	1
50-50-0300	CPR SEAT - FIXED	1
50-50-0510	SQUAD BENCH - BASE & LID WITH SHARPS AND TRASH	1
50-60-0002	OPTIONAL PATIENT COMPARTMENT ITEMS	1
50-60-0015	ASPIRATOR / SUCTION - IMPACT# 324 - SELF CONTAINED - COUNTER	1
50-60-0037	ASSIST HANDLES, (6) BLACK RUBBER	1
50-60-0153	COT FASTENER - STRYKER SINGLE - CENTER POSITION	1
50-60-0220	EMBLEMS / SIGNS - REAR PATIENT COMPARTMENT	1
50-60-0254	EXHAUST VENTILATOR, 3-SPEED - MULTIPLEX ELECTRIC SYSTEM	1
50-60-0300	FLOOR CONSTRUCTION- PATIENT COMPARTMENT, NON-WOOD COMPOSITE	1
50-60-0301	FLOORING - COVE	1
50-60-0340	GRAB RAIL, 1-1/4" DIAMETER, (96") OVER COT AREA	1
50-60-0390	HEADLINER - REAR PATIENT COMPARTMENT	1
50-60-0424	HEAT/AIR CONDITIONING, HOSELINE FT-454 4X4 AC/HEAT CORE	1
50-60-0434	I V HANGERS, DUAL CEILING, W/ VELCRO STRAPS - PERKO #1309	1
50-60-0555	LONPLATE II, NON-SLIP EMBOSSED SM. GRID - GUNMETAL #424	1
50-60-0612	OXYGEN OUTLETS (3) WALL - OHIO MEDICAL TYPE	1
50-60-0652	OXYGEN SYSTEM, MANUAL SYSTEM - MULTIPLEX ELECTRIC	1
50-60-0660	PADDED EDGING PROTECTION	1
50-60-0661	PADS, HEAD AND BACK	1
50-60-0700	SHARPS AND TRASH CONTAINERS - SQUAD BENCH END	1
50-60-0731	DOOR PANELS - MODULE PASSAGE DOORS - VINYL	1

