




## E-Z O<sub>2</sub> LIFT

Let's face it – oxygen bottles are heavy! The new E-Z O<sub>2</sub> Lift eliminates the heavy lifting by utilizing a nitrogen filled gas spring to raise and lower the bottle into the stowed position without electrical components. We have taken an age-old problem for many departments and created a solution that is cost effective, space saving and crew approved. To free up compartment space, the bottle is stored on the compartment door to maximize the storage area for life saving equipment.

### E-Z O<sub>2</sub> Lift Benefits:

- Smaller Footprint than Traditional Lift Systems
- Non Scissor Lift
- No Bottle Lifting
- Cost Savings vs. Other Lift Systems
- Reduce Back Injuries
- Less Downtime
- Ease of Use
- Non Electrical
- Reduces Compartment Damage when Stowing Bottles

The logo for the E-Z O<sub>2</sub> Lift. The letters "E-Z" are in a bold, red, sans-serif font. The "O<sub>2</sub>" is in a black, sans-serif font with a green "2" as a subscript. The word "Lift" is in a red, sans-serif font. The logo is positioned at the bottom of a large, detailed line drawing of the lift mechanism, which shows a green oxygen cylinder being raised and lowered by a white metal frame with a gas spring.

Scan for more info.