

QUOTATION

Braun Industries, Inc.

BRAUN INDUSTRIES, INC

CHIEF XL TYPE I 2019  
 FORD F450 193" WHEEL BASE 4X2  
 OVERALL L: 307" OVERALL W: 98" OVERALL H: 112"  
 MODULE L: 169" MODULE W: 96"  
 HEADROOM: 72"

Quote No: 10434-0001

AMB I: 19CHXLI

Chief XL 2019 Type I MOD 169"

02/13/2019

Page 1

PART NO	REF. NO	DESCRIPTION	QTY
		<b>== Chief XL 2019 Type I MOD 169" - 19.200 02/07/19 ==</b>	<b>1</b>
00-01-0006	STD CRASH TESTS	DYNAMIC SIDE IMPACT & ROLLOVER TESTS	1
	<b>CN8 COMPLIANCE REQ</b>	<b>KKK CHANGE NOTICE 8 REQUIRED</b>	<b>1</b>
00-01-0018	KKK-A-1822F CN8	BIDDER KKK-A-1822F CHANGE NOTICE 8 COMPLIANT	1
00-01-0019	KKK-A-1822F CN9	BIDDER KKK-A-1822F CHANGE NOTICE 9 COMPLIANT	1
00-01-0020	KKK-A-1822F CN10	BIDDER KKK-A-1822F CHANGE NOTICE 10 COMPLIANT	1
00-01-0021	STD WGT RATINGS CN10	-- WEIGHT RATINGS SPECIFIED, INTERIOR, LABELS, CN10 COMPLIANT	1
00-01-0022	KKK-A-1822F CN11	BIDDER KKK-A-1822F CHANGE NOTICE 11 COMPLIANT	1
00-01-0031		NO PERFORMANCE BOND REQUIRED	1
00-01-0041		NO BID BOND REQUIRED	1
00-01-0051		NO REQUIRED DELIVERY DATE	1
00-01-0053		PRECONSTRUCTION MEETING NOT REQUIRED / REQUESTED	1
00-01-0062	2017-XX LEAD TIME	LEAD TIME DEFINED	1
00-01-0091	STD 2015	REQUIRED REMAINING EQUIPMENT PAYLOAD:	1
	<b>SECTION HEADING</b>	<b>CHASSIS PACKAGES</b>	<b>1</b>
09-10-0001	OEM CHASSIS NOTE	<b>** IMPORTANT ** CHASSIS NOTIFICATION</b>	1
00-01-0223	2019 CHXL F450/550 4X2	-- TYPE I - AMBULANCE - CHIEF XL I - 72" INTERIOR HEADROOM	1
10-01-0194	2020 F450 193" WB 4X2	-- CHASSIS - FORD F-450 CHASSIS, 193" WB 4 x 2	1
10-08-0000	STD	-- ** NO CUSTOM CHASSIS REQUIRED **	1
30-10-0030	STD	-- BACKUP ALARM	1
30-10-0753	STD	-- SIREN, WHELEN 295-HFSA7-REMOTE WITH MICROPHONE	1
30-10-0812	STD 2017-XX F450	-- SIREN SPEAKERS, CAST PRODUCTS - IN FRONT BUMPER	1
20-10-0099	STD 2014	-- FUEL TANK - FULL UPON FACTORY RELEASE	1
20-10-0103	STD 2017 -XX CHFCLI	-- SUSPENSION - LIQUID SPRING	1

PART NO	REF. NO	DESCRIPTION	QTY
20-10-0206	STD F-SERIES	-- BATTERY - TWO SYSTEM - F SERIES CHASSIS	1
20-10-0603	STD 2014/15	-- ENGINE BLOCK HEATER - OEM PLUG	1
20-20-0134	STD 2019-XX F350/450	-- FLOOR CONSOLE - ANGLED FACE REMOVABLE PLATES - ALUMINUM	1
20-20-0252	STD 2018-XX	-- HEAT SHIELDS, PROTECTIVE	1
20-20-0277	STD E & F CHASSIS	-- HIGH IDLE SYSTEM, AUTOMATIC	1
20-30-0202	STD F SERIES	-- MIRRORS, EXTERIOR OEM MANUAL TELESCOPING W/REMOTE & HEAT	1
20-30-0261	STD F-SERIES	-- OEM AUTO-LOCK DISABLED (cab & module doors do not auto lock or unlock)	1
20-30-0806	STD 2014/15 TYPE I	-- RUNNING BOARDS - STAR PUNCHED DIAMOND PLATE	1
20-40-0206	STD	-- WHEEL SIMULATORS, STAINLESS STEEL WITH VALVE STEM EXTENSIONS	1
30-06-0350	STD 2014/15 M-TECH	-- FLASHING HEADLIGHTS - MULTIPLEX ELECTRICAL SYSTEM	1
30-25-3110	STD E-F-HD CHASSIS	-- LIGHTS, FRONT GRILLE - (2) WHELEN 500 SUPER LED - RED W/CLEAR LENS	1
30-27-1114	STD 2014/15 15D	-- LIGHTS, INTERSECTION FRONT - WHELEN 700 SUPER LED - RED W/CLR LENS	1
39-10-0003	STD TYPE I (12)	-- CONNECTIONS, MODULE	1
	<b>SECTION HEADING</b>	<b>ELECTRICAL SYSTEM FEATURES</b>	<b>1</b>
30-10-0050	STD	-- BATTERY CONVERTER / CHARGER	1
30-06-0004	STD 2017-XX	-- ELECTRICAL SYSTEM - WELDON V-MUX, MICRO PROCESSOR BASED	1
30-06-0010	STD 2014/15 M-TECH	-- VISTA IV MCC - PUSH BUTTON	1
30-06-0015	STD 2014/15 M-TECH	-- VISTA IV ACP - PUSH BUTTON	1
30-06-0026	STD 2014/15 M-TECH	-- WARNING DISPLAY - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0050	STD 2014/15 M-TECH	-- DAYTIME RUNNING LAMPS - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0076	STD 2014/15 M-TECH	-- MASTER CONTROL CONSOLE (MCC) - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0101	STD 2014/15 M-TECH	-- ATTENDANT CONTROL PANEL (ACP) - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0126	STD 2014/15 M-TECH	-- ELECTRICAL COMPARTMENT (PDQ) - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0150	STD 2014/15 M-TECH	-- SPARE CIRCUIT - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0175	STD 2014/15 M-TECH	-- WIRING HARNESS CONNECTORS - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0201	STD 2018-XX	-- RADIO TIE-IN POINT - (1) 12V/40 AMP - CAB CONSOLE AREA	1
30-06-0202	STD 2018-XX	-- RADIO TIE-IN POINT - (1) 12V/40 AMP - PATIENT COMPARTMENT	1
30-06-0203	STD 2018-XX	-- RADIO TIE-IN POINT - (1) 12V/40 AMP - PDQ	1
30-06-0225	STD 2014/15 M-TECH	-- EMERGENCY MASTER SWITCH - MULTIPLEX ELECTRICAL SYSTEM	1

PART NO	REF. NO	DESCRIPTION	QTY
30-06-0250	STD 2014/15 M-TECH	-- CONTROL PANEL SWITCHES - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0275	STD 2014/15 M-TECH	-- MODULE POWER SWITCH - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0285	STD 2018-XX M-TECH	-- MODULE DISCONNECT - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0300	STD 2014/15 M-TECH	-- VOLTMETER - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0325	STD 2014/15 M-TECH	-- FLASHER - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0375	STD 2014/15 M-TECH	-- PARK OVERRIDE - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0400	STD 2014/15 M-TECH	-- LOAD MANAGER-SEQUENTIAL SWITCHING SYSTEM - MULTIPLEX ELECTRIC SYSTEM	1
30-06-0410	STD 2016-XX M-TECH	-- WARNING LIGHTS FLASH PATTERN SIGNAL ALERT 75 - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0450	STD 2014/15 M-TECH	-- FRONT CLEAR DISABLE SWITCH - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0475	STD 2014/15 M-TECH	-- CURBSIDE SCENE LIGHTS "ON" WITH CURBSIDE DOOR - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0500	STD 2014/15 M-TECH	-- REAR SCENE LIGHTS "ON" WITH REAR DOORS OPEN - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0525	STD 2014/15 M-TECH	-- DISABLE SWITCH FOR REAR SCENE LIGHTS - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0550	STD 2014/15 M-TECH	-- DISABLE SWITCH FOR SIDE SCENE LIGHTS - MULTIPLEX ELECTRICAL SYSTEM	1
30-06-0575	STD 2014/15 M-TECH	-- REVERSE (TRANSMISSION) FUNCTIONS - MULTIPLEX ELECTRICAL SYSTEM	1
30-10-0574	STD 2019-XX	-- USB - DUAL PORT OUTLETS, TWO (2)	1
30-10-0901	STD WITH V-MUX	-- SYSTEM ACTIVE AND CHECK OUT LIGHTS - 15 MINUTES	1
30-11-0106	STD 2018-XX LED DOME	-- LIGHTS, DOME SUPER LED (8) WHELEN	1
50-60-0256	STD 2014/15	-- EXHAUST VENTILATOR, 3-SPEED - MULTIPLEX ELECTRIC SYSTEM	1
50-60-0640	STD 2014/15	-- OXYGEN SYSTEM, ELECTRIC - MULTIPLEX ELECTRIC SYSTEM	1
50-60-0629	STD	-- FLOW METER, OXYGEN - THORPE STYLE	1
30-10-0071	STD 2013 CHXL/SCHF	CAMERA REAR VIEW - SINGLE SYSTEM - VISTA SCREEN	1
30-10-0110	STD	ELECTRICAL CIRCUITS	1
30-10-0190	STD 2014/15	GROUND FAULT CIRCUIT INTERRUPTER (GFCI)	1
30-10-0291	STD 2017-XX CXL/SCHF	LIGHT, ATTENDANT CONTROL PANEL AREA - LED - WARM WHITE	1
30-10-0310	STD CHXL/SCHF/TLC	LIGHT, LED ELECTRICAL COMPARTMENT (PDQ)	1
30-10-0410	STD CHXL /SCH 10 LTS	LIGHTS, OSS COMPARTMENTS - INCANDESCENT	1
30-10-0527	STD	OUTLETS, 125VAC - THREE (3) INTERIOR DUPLEX	1
30-10-0543	STD	OUTLETS, 12VDC - TWO (2) INTERIOR POWERPOINT	1
30-10-0657	STD TYPE I	RADIO ANTENNA BASES - TWO (2)	1
30-10-0723	STD -NO BLOCK HEA	RECEPTACLE SHORELINE INLET-SUPER 20 AMP, AUTO EJECT	1
30-11-0152	STD 2019-XX	LIGHTS, LED PANEL (3) STREETSIDE - ANGLED	1
30-11-0176	STD 2019-XX	LIGHTS, LED PANEL (2) CURBSIDE - ANGLED	1
<b>SECTION HEADING</b>		<b>EXTERIOR LIGHTING</b>	<b>1</b>

PART NO	REF. NO	DESCRIPTION	QTY
30-25-0122	STD - TECNIQ	-- LIGHTS, ICC / LED CLEARANCE - ON MODULE FRONT - TECNIQ	1
30-25-4201	STD 2014/15 CHXL/SCHF	-- LIGHTS, FRONT MODULE - WHELEN (7)900 - (6)RED, (1) CLEAR LED W/CLEAR LENS	1
30-26-3012	STD 2019-XX	LIGHTS, REAR TAIL - WHELEN 600 LED STOP, TURN, & MINIMUM INTENSITY BACK-UP	1
30-26-5101	STD 2014/15 CHXL/SCF/S	LIGHTS, REAR WINDOW-LEVEL - WHELEN (2)900 LED - RED W/CLR LENS	1
30-25-0123	STD SS - TECNIQ	-- LIGHTS, ICC / LED CLEARANCE - ON MODULE REAR - TECNIQ	1
30-26-6155	STD 2018-XX	-- LIGHTS LED, REAR MODULE -WHELEN (5)900: (2)RED, (1)AMB , (2)SCN CLR LENS	1
30-26-7002	STD 2009	LIGHT, LICENSE PLATE LED	1
30-27-1102	STD - TECNIQ	LIGHTS, ICC / LED CLEARANCE - ON MODULE SIDE - TECNIQ	1
30-27-2004	STD 2018-XX	LIGHTS, SIDE MODULE SCENE - WHELEN 900 24 DIODE LED - CLEAR	1
30-27-3009	STD 2018-XX C/S/T	LIGHTS, SIDE MODULE WARNING - WHELEN 900 SUPER LED - RED W/CLR LENS	1
30-27-4115	STD 2014/15 CHXL SCHF	LIGHTS, INTERSECTION REAR - WHELEN 700 SUPER LED - RED W/CLR LENS	1
	<b>SECTION HEADING</b>	<b>MODULE BODY CONSTRUCTION AND CONNECTIONS</b>	<b>1</b>
39-00-0001	STD SOLIDBODY	MODULE BODY CONSTRUCTION - SEAMLESS	1
39-00-0002	STD BODY CONSTRUCTION	SOLE SOURCE MODULE BODY CONSTRUCTION	1
39-00-0005	STD ALL	AWS ALUMINUM CERTIFIED WELDERS	1
	<b>SECTION HEADING</b>	<b>MODULE FEATURES / EXTERIOR</b>	<b>1</b>
40-00-0009	STD 2014	ALUMINUM PARTS - MODULE EXTERIOR	1
40-00-0010	STD	CORROSION MINIMIZATION - MODULE	1
40-00-0019	STD	DOOR SKINS, OUTSIDE STORAGE COMPARTMENTS - DIAMOND PLATE INTERIOR SKINS	1
40-00-0020	STD	DOOR CONSTRUCTION - GASKET PLACEMENT	1
40-00-0030	STD	EXTERIOR COMPARTMENTS DESIGN	1

PART NO	REF. NO	DESCRIPTION	QTY
40-00-0050	STD	EXTERIOR COMPARTMENTS, VENTING	1
40-00-0055	STD	EXTERIOR COMPARTMENTS, VENTING - OXYGEN STORAGE	1
40-00-0060	STD	EXTERIOR COMPARTMENTS - SWEEP-OUT	1
40-00-0070	STD	HANDLES, EBERHARD - CHROME FREE FLOATING TYPE	1
40-00-0080	STD	INSULATION - POLYURETHANE SPRAY FOAM	1
40-00-0091	STD TYPE I	THRESHOLD, STAINLESS STEEL - ALL OSS STORAGE COMPARTMENTS	1
40-10-0066	STD 2014/15	-- DOOR SLIDING - MODULE CURBSIDE PASSAGE	1
50-60-0037	STD EXP/SS/RESP	-- ASSIST HANDLES, (6) BLACK RUBBER	1
40-10-0072	STD 2015	LOWERED MODULE SIDES	1
40-10-0088	STD 2018-XX	DOORS - REAR PASSAGE	1
40-10-0103	STD 2014/15	FENDERETTES - POLISHED ALUMINUM	1
40-10-0120	STD 2014/15	HOLD-OPEN - REAR DOORS - GRABBERS	1
40-10-0130	STD 2018-XX ALL	EMERGENCY RELEASE LEVER, EXPOSED, REAR PATIENT PASSAGE DOORS	1
40-10-0189	STD2014/15 5 OSS-OEM F	LOCKS POWER - PASSAGE, COMPARTMENTS, KEY FOBS	1
40-10-0151	STD 2017-XX UNLOCK	-- SWITCH UNLOCK, WATERPROOF - HIDDEN IN FRONT GRILLE AREA	1
40-10-0204	STD 2015	MUD FLAPS, REAR - BLACK RUBBER FIBERGLASS REINFORCED - THREE COLOR	1
40-10-0266	STD 2014/15 STAR PUNCH	REAR STEP/BUMPER -FLIP-UP CENTER SECTION - STAR PUNCHED AND ANODIZED	1
40-10-0280	STD 2014	RUB RAILS, ANODIZED ALUMINUM	1
60-25-0101	STD	-- BLACK REFLECTIVE MATERIAL IN RUB RAILS	1
40-10-0321	STD 2014/15	STEPWELL, CURBSIDE DOORWAY - MID STEP AREA	1
40-10-0340	STD 2014/15	STONE GUARDS, MODULE FRONT - ANODIZED DIAMOND PLATE	1
40-10-0360	STD 2014/15	KICKPLATE/THRESHOLD, REAR DOORWAY - ANODIZED DIAMOND PLATE	1
40-10-0381	STD 2014/15 CHXL/LIB	TOW HOOKS, MODULE REAR	1
40-10-0505	STD	WINDOWS, PASSAGE DOOR - (CURB SLIDING, REAR FIXED)	1
42-01-0202	STD 2017-XX CHXL I D-F	OUTSIDE STORAGE #2 - STREETSIDE CENTER	1
42-02-0107	STD 2017-XX	-- SHELF ONE (1) IN OSS #2	1
42-02-0056	STD 2017-XX OSS #2	-- SHELF TRACK - EXPOSED - OSS #2	1
43-01-0204	STD 2015 CHXL I	OUTSIDE STORAGE #3 - STREETSIDE REAR - TALL DOOR	1
43-02-0108	STD 2015	-- SHELF ONE (1) IN OSS #3	1
43-02-0051	STD 2015	-- SHELF TRACK - EXPOSED - OSS #3	1
41-01-0207	STD 2015 CHXL I D-F	-- OUTSIDE STORAGE #1 - STREETSIDE FORWARD	1
44-01-0208	STD 2018-XX CHXL I	-- OUTSIDE STORAGE #4A BOARD/SCOOP, #4B O2 CYLINDER - RT. REAR	1
50-30-0114	STD 2018-XX CHFCLI	-- CABINET - RIGHT WALL, REAR WITH (1) DOOR FACING AISLE #1	1
45-01-0203	STD 2015 CHXL I	OUTSIDE STORAGE #5 - CURBSIDE FORWARD	1
<b>SECTION HEADING</b>		<b>MODULE INTERIOR / PATIENT COMPARTMENT FEATURES</b>	<b>1</b>

PART NO	REF. NO	DESCRIPTION	QTY
50-00-0202	STD 2014/15 CH/SCH/	ROUNDED INTERIOR CORNERS	1
50-51-1001	STD	-- CABINET CONSTRUCTION - INTERIOR - ALUMINUM	1
50-51-1011	STD ALUM CABINETS	-- RECESSED ADJUSTABLE SHELF TRACK - CABINETS	1
50-51-1013	STD 2017-XX	-- HOLD OPENS FOR INTERIOR DOORS	1
50-51-1014	STD 2016-XX	-- CATCH, MAGNETIC - ADHESIVE BACK FOR INTERIOR DOOR(S) OVER 36 INCHES	1
50-51-2002	STD 2018-XX	-- CABINET DOORS - CLEAR POLYCARBONATE	1
50-51-3001	STD LB/CH/SC/TLC	-- CABINETS/ WALLS - POLYCHROMATIC COATING - MARBLE STONE	1
50-51-4012	STD 2018 CHXL/SCHF	-- COUNTER TOP - SOLID ACRYLIC - NIGHT STARS (BLACK)	1
50-51-5001	STD 2018-XX	-- UPHOLSTERY - PATRIOT PLUS CHARCOAL #8605	1
50-51-6003	STD CHXL/BDD/SCF	-- LONPLATE II, NON-SLIP EMBOSSED SM. GRID - GUNMETAL #424	1
50-60-0391	STD 1 PC ALUM WELDED	-- HEADLINER - REAR PATIENT COMPARTMENT	1
50-01-0408	2018-XX LATCHES	LATCHES, CABINET, CN10 COMPLIANT AS APPLICABLE	1
50-10-0110	2018 STD CHF XL	-- CABINET - LEFT WALL, UPPER REAR #1 - FIXED FRAME	1
50-10-0229	STD 2016-20XX	-- LEFT WALL - ATTENDANT CONTROL PANEL - HINGED PANEL	1
50-10-0237	STD 2018 CHF XL	-- CABINET - LEFT WALL, UPPER OVER PANEL #2 - FIXED FRAME	1
50-10-0340	STD 2018 CHXL I	-- CABINET - LEFT WALL, LOWER REAR #3 - INSIDE/OUTSIDE	1
50-10-0430	STD 2018 XX CHXL	-- CABINET - LEFT WALL, REAR ACTION AREA #4	1
50-10-0604	STD 2018 CHXLI & BDDI	-- PRIMARY ACTION AREA	1
50-10-0637	2018 - 20XX STD	-- SECONDARY ACTION AREA CHF XLI	1
50-10-0705	STD 2014/15	-- DRAWER - SLIDE-OUT / WRITING TABLE	1
50-50-0324	STD 2019-20XX CHF/SCH	STREETSIDE INTERIOR  -- CPR SEAT FIXED, VALOR A-VAL BACK, CN8 COMPLIANT	1
50-20-0147	STD 2018-XX LIB I	-- CABINET - FRONT WALL, UPPER #1	1
50-20-0289	STD 2018-XX CHF XLI	-- CABINET - FRONT WALL, LOWER #2 & #3 - INSIDE/OUTSIDE ACCESS	1
50-20-0381	STD CHXL I	-- CABINET - FRONT WALL, LOWER BELOW PDQ (SINGLE VENTED DOOR)	1
50-20-0507	STD	-- PASS THROUGH - CAB TO MODULE, WINDOW OPENING	1
50-41-0213	STD 2019- XX	-- HEAT/AC PROAIR IW-3856 580 CFM 4X4 HEAT/AC CORE, BRUSHLESS BLOWER MOTOR	1

PART NO	REF. NO	DESCRIPTION	QTY
50-50-0118	2018-XX STD CN8	ATTN SEAT - USSC VALOR R - BACK INTEGRAL CHILD, SWIVEL SEAT CN8 COMPLIANT	1
50-50-0028	STD 2019-XX CHF/SCHF	-- SEATING - SQUAD BENCH BASE W/ STORAGE CN8 COMPLIANT VALOR A-VAL BACK	1
50-50-0503	STD 2018-XX	-- SQUAD BENCH TOP - SPLIT WITH HINGED LID FORWARD WITH SHARPS / TRASH ACCESS	1
50-50-0541	STD	-- HOLD OPENS - SQUAD BENCH LID - GAS CYLINDERS	1
50-60-0709	STD 2014/15	-- SHARPS AND TRASH CONTAINERS - ACCESS THRU SQUAD BENCH LID	1
50-60-0718	STD 2014/15	-- SQUAD BENCH FACE - VINYL FLOORING MATERIAL	1
	<b>2018 CHFXL COT MNT FAS</b>	<b>COT MOUNT FASTENERS</b>	<b>1</b>
50-60-0161	STD 2018-XX 96"	-- COT FASTENER - STAT-TRAC SINGLE - CENTER POSITION CN8 COMPLIANT	1
	<b>2014/15 CHXL I</b>	<b>OPTIONAL PATIENT COMPARTMENT ITEMS</b>	<b>1</b>
50-60-0018	STD RECESSED 2016 -XX	ASPIRATOR / SUCTION - SSCOR - WITH ELECTRIC PUMP - RECESSED	1
50-60-0220	STD	EMBLEMS / SIGNS - REAR PATIENT COMPARTMENT	1
50-60-0300	STD	FLOOR CONSTRUCTION- PATIENT COMPARTMENT, NON-WOOD COMPOSITE	1
50-60-0301	STD	FLOORING - COVE	1
50-60-0322	STD CHXL/SCF/PAT	GLOVE DISPENSER - (3) BOX - ABOVE ENTRY DOOR - CURBSIDE (ABS)	1
50-60-0340	STD	GRAB RAIL, 1-1/4" DIAMETER, (96") OVER COT AREA	1
50-60-0431	STD 2018-XX	I V HANGERS - CEILING RECESSED - CAST PRODUCTS #IV2008-1 (2)	1
50-60-0613	STD SS/CXL/SCHF	OXYGEN OUTLETS (2) WALL/(1)CEILING - OHIO MEDICAL TYPE	1
50-60-0660	STD	PADDED EDGING PROTECTION	1
50-60-0661	STD	PADS, HEAD AND BACK	1
50-60-0730	STD 2014/15	DOOR PANELS - MODULE PASSAGE DOORS - FULL LENGTH ALUMINUM - W/ACCESS PANELS	1
50-60-0761	STD 2014/15	STAINLESS STEEL WALL PROTECTION - INTERIOR STREETSIDE	1
50-60-0820	STD 2014/15	TURTLE TILE - SKID-RESISTANT MAT, CURBSIDE STEPWELL	1
	<b>SECTION HEADING</b>	<b>PAINT / GRAPHICS</b>	<b>1</b>
60-01-0001	STD	PAINT PROCEDURE	1
60-01-0002	STD	PAINT MODULE - OEM WHITE	1
60-01-0003	STD	CLEAR COATING - MODULE PAINT	1
60-01-0004	STD	BUFFING, PAINT - MODULE	1

PART NO	REF. NO	REF. NO	DESCRIPTION	QTY
60-10-0108	STD		CAB - OEM PAINT	1
60-30-0200	STD 2017-XX		DOOR REFLECTIVITY	1
80-10-0000	STD		KKK-A-1822F COMPLIANCE	1
	<b>STD</b>		<b>VEHICLE MANUALS</b>	1
80-20-0001	STD		VEHICLE MANUALS (1) PACKAGE SET	1
90-10-0007	2018 STD WARRANTY NEW		STATEMENT OF LIMITED WARRANTY NEW AMBULANCES	1