



AMBULANCE REMOUNT PACKAGES

Maximize your fleet's value with our three comprehensive ambulance remount solutions. Each package is customizable to meet the specific needs of your department or organization and has the opportunity to extend your factory warranty (see website for more details). Contact your local Braun ambulance dealer to help select the right package for you.

1. REACTIVATE REMOUNT PACKAGE

For departments that need a fast and focused return to service.

- Best suited for departments seeking a quick chassis replacement
- Includes limited add-on options
- Module placed on a new chassis in the same make and model (or a very close match) to the original chassis model



2. REVIVE REMOUNT PACKAGE

Breathe new life into your existing module with core updates and refinishing.

- Designed for aging ambulances needing updates with limited customization options available
- Module placed on a new chassis in the same make and model (or a very close match) to the original
- Includes anodized aluminum diamond plate replacement, module repaint, upholstery repair, and floor replacement



3. REFRESH REMOUNT PACKAGE

Transform your current module into a modern, fully equipped frontline unit.

- Our most comprehensive remount solution with the greatest flexibility
- Includes all upgrades from package #2, plus expanded options
- Ideal for departments wanting a like-new ambulance



For the protection of the remount team, the module must be cleaned and sanitized (including the removal of waste & sharps) prior to customer/department drop off at the remount facility.

AMBULANCE REMOUNTS

SERVICES & FEATURES INCLUDED BY PACKAGE:

| | 1 | 2 | 3 |
|---|---|---|---|
| New Body Mounts | ✓ | ✓ | ✓ |
| New Cab Running Boards | ✓ | ✓ | ✓ |
| New Siren / Horns / Speaker | ✓ | ✓ | ✓ |
| Reinstall Screens & Replace Buttons | ✓ | ✓ | ✓ |
| HVAC Hook Up | ✓ | ✓ | ✓ |
| New HVAC Filter | ✓ | ✓ | ✓ |
| New Back Up Alarm | ✓ | ✓ | ✓ |
| New Grille Lights | ✓ | ✓ | ✓ |
| Transfer Radio Antenna Base Transfer | ✓ | ✓ | ✓ |
| Replace Grit Tape | ✓ | ✓ | ✓ |
| PDI / Clean Unit | ✓ | ✓ | ✓ |
| Repainted Module (White Standard - Charge for Optional Color) | | ✓ | ✓ |
| New Rub Rails | | ✓ | ✓ |
| New Fuel Fill / DEF Fill | | ✓ | ✓ |
| New Rear Bumper | | ✓ | ✓ |
| New Diamond Plate Protective Panels | | ✓ | ✓ |
| New Fender Flares | | ✓ | ✓ |
| New Drip Rail | | ✓ | ✓ |
| New Vinyl Floor | | ✓ | ✓ |
| New Module Switch Panel & Switches | | ✓ | ✓ |
| Service Compartment / Entry Door Latches | | ✓ | ✓ |
| Repair Upholstery if Possible (ID at Inspection) | | ✓ | ✓ |
| Limited Optional Add-ons | ✓ | ✓ | |
| Full Optional Add-ons | | | ✓ |

OPTIONAL ADD-ONS BY PACKAGE:

| | 1 | 2 | 3 |
|--|---|---|---|
| LiquidSpring Suspension | ✓ | ✓ | ✓ |
| Power Door Locks / New Stealth Switch | ✓ | ✓ | ✓ |
| Transfer Communication Equipment (Final Radio Hook-ups Not Included) | ✓ | ✓ | ✓ |
| Graphics | ✓ | ✓ | ✓ |
| Paint Cab to Match Custom Module Paint | ✓ | ✓ | ✓ |
| Replace Plexi - Same Latches | | | ✓ |
| Replace Upholstery - Subject to Availability | | | ✓ |
| Re-Spray Exterior Compartment Tubs | | | ✓ |
| Update HVAC (Match Original System) | | | ✓ |
| Upgrade Cot Fastener | | | ✓ |
| Upgrade Exterior Lights to LED | | | ✓ |
| Upgrade Dome Lights to LED | | | ✓ |
| Replace Exterior Compartment/Entry Door Handles (With Same Model Received) | | | ✓ |
| Replace Interior Door Handles/Latches | | | ✓ |



START YOUR AMBULANCE REMOUNT PROJECT AT
BraunAmbulances.com



1170 Production Drive
 Van Wert, OH 45891

Toll Free (877) 344-9990